MUX8-R2™

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Compatible with
PalmSens3, (Multi)PalmSens4, EmStat3 Blue and
EmStat4X LR and HR

➤ See for more information:
www.palmsens.net

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
MUX8-R2 Settings

- Connect Sense to WE
- Combine RE and CE
- Use Common RE and CE on Channel 1

Unselected WE

- Disconnect WE (floating)
- Switch WE to GND

The hardware configuration can be set in our software.

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
Cell Connection Options

Option A (default):

The channels are divided in two sets of four sensor cables joined with a D-Sub connector.

Order code: [CBL-MUX08R2-SNS-5S]

Option B:

The cable here shown at the right can be used in case the multiplexer needs to be connected to a fixed setup by means of soldering or screw-terminals.

Order code: [CBL-MUX08R2-HD]

Option C:

You can also connect one or two screw-terminals directly in the multiplexer.

Order code: [MUX08R2-ST]

Option D:

The SPE adapter for our MUX8-R2 multiplexer allows you to connect 8x Screen Printed Electrodes (SPE's). The pitch of the SPE connector is 2.54 mm and compatible with the most popular brands of SPE's.

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Voltammetric techniques	Sequential	Alternating
▪ Linear Sweep Voltammetry	✓	-
▪ Cyclic Voltammetry	✓	-
▪ Fast Cyclic Voltammetry	✓	-
▪ AC Voltammetry	✓	-
▪ Differential Pulse Voltammetry	✓	-
Pulsed techniques:		
▪ Square Wave Voltammetry	✓	-
▪ Normal Pulse Voltammetry	✓	-
▪ Stripping Chronopotentiometry	✓	-
Amperometric techniques		
▪ Chronoamperometry	✓	✓
▪ Zero Resistance Amperometry	✓	✓
▪ Multistep Amperometry	✓	-
▪ Fast Amperometry	✓	-
▪ Pulsed Amperometric Detection	✓	-
▪ Multiple-Pulse Amperometric Detection	✓	-
Galvanostatic techniques		
▪ Linear Sweep Potentiometry	✓	-
▪ Chronopotentiometry	✓	✓
▪ Multistep Potentiometry	✓	-
▪ Open Circuit Potentiometry	✓	✓
▪ Stripping Chronopotentiometry	✓	-
Other		
▪ Mixed Mode	✓	-
▪ Impedance Spectroscopy (EIS/GEIS)	✓	-

Specifications


System specifications

▪ number of channels	8 (up to 128 channels when daisy chained)
▪ multiplexer	switches 8 x [WE, WE2 or Sense, RE, CE]
▪ on resistance for WE	1.5 ohm typical
▪ charge injection for WE	20 pC typical
▪ leakage current	< 55 pA (20 pA typical) at 25 °C
▪ switching time	2 ms
▪ compliance voltage	±10 V
▪ maximum current	70 mA (when used with EmStat3+ Blue or EmStat4X HR)


Limitations for Electrochemical Impedance Spectroscopy (EIS)

▪ max. frequency for EIS	100 kHz when switching WE/S, RE and CE 1 MHz when switching WE/S and RE+CE combined (2 electrodes configuration)
--------------------------	--

Functional Diagram

Dimensions

Dimensions in mm

Software Support

The MUX8-R2 is supported by PStouch and MultiTrace for Windows and the PStouch app for Android.

MUX8-R2 with Integrated Potentiostat

EmStat MUX8-R2™

The MUX8-R2 multiplexer is also available with integrated EmStat3 or EmStat4 potentiostat. This very compact combination allows for a high productivity with a small footprint.

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DISCLAIMER

Changes in specifications and typing errors reserved. Every effort has been made to ensure the accuracy of this document. However, no rights can be claimed by the contents of this document.